

Interim PALS CPR/AED Testing Checklist Child 1-Rescuer CPR and AED Test

Released: November 11, 2010

Student Name:	Date of Test:	
---------------	---------------	--

Skill Step	Critical Performance Steps	Child CPR With AED ☑ if done correctly	
1	Checks unresponsiveness		
2	Checks for no breathing or only gasping		
	Minimum 5 seconds; maximum 10 seconds		
3	Activates emergency response system/calls for AED	Y	
4	Checks carotid pulse		
	Minimum 5 seconds; maximum 10 seconds		
5	Locates CPR hand position		
6	Delivers first cycle of compressions at correct rate		
	Acceptable <18 seconds for 30 compressions		
7	Gives 2 breaths (1 second each)		
8	Delivers second cycle of compressions at correct hand position		
	Acceptable >23 out of 30 compressions		
9	Gives 2 breaths (1 second each)		
	AED arrives		
AED	Turns AED on		
1			
AED	Selects proper AED pads and places correctly		
2			
AED	Clears victim to analyze		
3	(Must be visible and verbal check)		
AED	Clears victim to shock/presses shock button		
4	(Must be visible and verbal check)		
Ctudon	Maximum time from AED arrival 90 seconds t continues CPR		
Studen 10/	Delivers third cycle of compressions at correct rate		
10	Acceptable <18 seconds for 30 compressions		
11	Gives 2 breaths (1 second each) with visible chest rise		
	t step is done only with a manikin with a feedback device, such as a click	er or light If no	
	k device, STOP THE TEST.	er or ugia. If no	
12	Delivers fourth cycle of compressions of adequate depth with full		
	chest recoil		
	Acceptable >23 compressions		
	•	STOP THE TEST	
		PASS NR	