

Client Memo

To: Wake Forest Baptist Health Outreach Clients
Date: July 7, 2017
Re: New Tube for Troponin Testing

On July 7 2017, the Clinical Laboratory will begin using a new tube type for Troponin levels. The test will still appear in CareEvolve as **Troponin** (test code **LAB747**). The required tube type is the only thing that is changing..

- The new container code will be: Green Top-Li/H.
- The container description will be: GrnLi/H.
- The minimum volume will be: 25 mL.

The Clinical Laboratory is switching to the new specimen container test for Troponin levels. Although, a gold top or red top will still be accepted until Thursday, July 13, 2017, after that the only acceptable container will be the Green Top Li/H.

If you have any questions related to this test change, please contact the Core Laboratory at 716-2610.